

Scanning Tunneling Spectroscopy (STS) experiments were performed to demonstrate the ultimate performance of our CreaTec LT-STM systems. Therefore, a superconducting Pb(111) single-crystal and a Pb tip were used to measure STS spectra at 1.8 K using the latest CreaTec DSP electronics. The modulation amplitude of the internal CreaTec Lock-In Amplifier was set to very small values in the range of a few  $\mu\text{V}$ , i.e. far below  $kT$ , to eliminate the influence of the modulation amplitude. As a consequence, the main contribution to the energy resolution is the temperature of the STM itself.

## SPECTRA




Fig. 1: I(V) spectrum clearly showing the superconducting gap in the tunneling current.




Fig. 2:  $dI/dV$  spectrum of the same sample at a different modulation amplitude.




Fig. 3:  $dI/dV$  spectra of the same sample at different modulation amplitudes in the range of 5-200  $\mu\text{V}$  (bottom to top).

The spectra are non-averaged single spectra (raw data) recorded with the following parameters:

- 1024 points/spectrum
- 15 second/spectrum
- Bias Voltage 10 mV
- Tunneling Current setpoint 100 pA
- $T=1.8 \text{ K}$  ( $kT=154 \mu\text{eV}$ )