

CreaTec offers custom designed water cooling systems for evaporators and MBE systems. CreaTec's water cooling systems are individually adapted to the effusion cell type and the requirements of the customer's MBE system.

Options

- Closed cooling shroud (mounted directly to the flange)
- Open cooling shroud (stands free, only the tubes are connected to the flange)
- Flow controller

TYPE	WK
MATERIAL	stainless steel
COOLING MEDIUM	water alcohol
FLOW RATE	2 - 5 l/min
PRESSURE	1 - 2 bar
FLANGE SIZE	DN 40 CF 63 100
DESIGN AND GEOMETRY	as specified
WATER CONNECTION	standard tube 6 x 1 or as specified

